

Exclusive for you

ROLLATOR

Find us:



X Folding - Simple handle - Light weight

The Rollators will be your faithful companion when you go outdoor. It's guarantee a very comfortable ride and its perfect for every urban and outdoor area.

Introducing our new range of centre folding rollators:

Rollator
RG78



NEW

Bariatric Rollator
RA66 BR



NEW

Narrow Rollator
RA60 N



NEW

Specification

Handle height	87.5 cm - 99.5 cm
Seat height	61 cm
Length	65 cm
Width	65 cm
Width between handle	48 cm
Net weight	6.3 kg
Weight capacity	136 kg / 300 lbs

Specification

Handle height	74 cm - 94 cm
Seat height	54 cm
Length	65 cm
Width	70 cm
Width between handle	53 cm
Net weight	6.5 kg
Weight capacity	205 kg / 450 lbs

Specification

Handle height	70 cm - 90 cm
Seat height	48 cm
Length	65 cm
Width	58 cm
Width between handle	43 cm
Net weight	7.0 kg
Weight capacity	136 kg / 300 lbs

Features:



Double crossbar



Folds easy for storage



Airmesh Fabric Seat

Features:



Large backrest



Easily Height adjustable handles



Airmesh fabric seat

Features:



Slow down brake



Folds easy for storage



Small roller designed

Trolli Master Sdn. Bhd. (209382-P)
Lot 525, Mukim Telok Panglima Garang,
42500, Daerah Kuala Langat
Selangor Darul Ehsan, Malaysia

Tel: +603 3191 4077 / +603 3191 2837
Fax: +603 3191 2152 / +6033191 3022
Email: tm@trollimaster.com.my
Website: www.trollimaster.com





Specification

	STANDARD	LOW	SUPER LOW
Handle height	81 - 101 cm	74 - 94 cm	70 - 90 cm
Seat height	61 cm	54 cm	48 cm
Length	65 cm	65 cm	65 cm
Width	65 cm	65 cm	65 cm
Width between handle	48 cm	48 cm	48 cm
Net weight	7.2 kg	7.2 kg	7.2 kg
Weight capacity	136 kg/300lbs	136 kg/300lbs	136 kg/300lbs

RA60

- 7 Handle height adjustments.
- Available with basket or shopping bag.
- Maximum basket and shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.



Specification

	STANDARD	LOW	SUPER LOW
Handle height	81 - 101 cm	74 - 94 cm	70 - 90 cm
Seat height	61 cm	54 cm	48 cm
Length	65 cm	65 cm	65 cm
Width	65 cm	65 cm	65 cm
Width between handle	48 cm	48 cm	48 cm
Net weight	6.2 kg	6.1 kg	6.0 kg
Weight capacity	136 kg/300lbs	136 kg/300lbs	136 kg/300lbs

RA66

- Cable hidden braking system
- 7 Handle height adjustments.
- Available shopping bag only.
- Maximum shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.



Specification

	STANDARD	LOW	SUPER LOW
Handle height	80 - 100 cm	72 - 92 cm	69.5 - 88.5 cm
Seat height	58 cm	51 cm	46 cm
Length	65 cm	65 cm	65 cm
Width	65 cm	65 cm	65 cm
Width between handle	48 cm	48 cm	48 cm
Net weight	6.9 kg	6.8 kg	6.7 kg
Weight capacity	125 kg/275lbs	125 kg/275lbs	125 kg/275lbs

RA62

- 7 Handle height adjustments.
- Available with basket or shopping bag.
- Maximum basket and shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.



Specification

	STANDARD	LOW	SUPER LOW
Handle height	87.5 - 99.5 cm	79.5 - 88.5 cm	76.5 - 85.5 cm
Seat height	61 cm	54 cm	48 cm
Length	65 cm	65 cm	65 cm
Width	65 cm	65 cm	65 cm
Width between handle	48 cm	48 cm	48 cm
Net weight	7.3 kg	7.2 kg	7.1 kg
Weight capacity	136 kg/300lbs	136 kg/300lbs	136 kg/300lbs

RG70

- 4 Handle height adjustments.
- Available with basket or shopping bag.
- Maximum basket and shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.



Specification

	STANDARD	LOW	SUPER LOW
Handle height	85.5 - 94.5 cm	77.5 - 86.5 cm	74 - 83 cm
Seat height	61 cm	54 cm	48 cm
Length	65 cm	65 cm	65 cm
Width	65 cm	65 cm	65 cm
Width between handle	48 cm	48 cm	48 cm
Net weight	6.9 kg	6.8 kg	6.7 kg
Weight capacity	136 kg/300lbs	136 kg/300lbs	136 kg/300lbs

RG72

- 4 Handle height adjustments.
- Available with basket or shopping bag.
- Maximum basket and shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.



Specification

	STANDARD	LOW	SUPER LOW
Handle height	87.5 - 99.5 cm	79.5 - 88.5 cm	76.5 - 85.5 cm
Seat height	61 cm	54 cm	48 cm
Length	65 cm	65 cm	65 cm
Width	65 cm	65 cm	65 cm
Width between handle	48 cm	48 cm	48 cm
Net weight	6.3 kg	6.2 kg	6.1 kg
Weight capacity	136 kg/300lbs	136 kg/300lbs	136 kg/300lbs

RG77

- 4 Handle height adjustments.
- Available shopping bag only.
- Maximum shopping bag load 5kg / 11lbs.
- 200mm EVA wheel.